


## HSV GFP Control Virus

*HindIII* *KpnI*  
 5'-- AAG CTT GGT ACC GAG CTC GGA TCC ACT AGT AAC

*EcoRI* *EcoRV*  
 GGC CGC CAG TGT GCT GGA ATT CTG CAG ATA TCC ATC

*NatI* *XhoI* *XbaI*  
 ACA CTG GCG GCC GCT CGA GCA TGC ATC TAG A --3'



Cat.	Titer	Quantity
C172	HSV GFP Control Virus (10 <sup>6</sup> IU/ml)	10 ml
C179	HSV GFP Control Virus (10 <sup>7</sup> IU/ml)	2 x 200 µl

<b>Species</b>	N/A
<b>Accession Number</b>	N/A
<b>Vector</b>	pHSV-CMV-GFP-Neo Control Vector
<b>Vector Size</b>	6055bp
<b>Insert Size</b>	720bp
<b>Sequencing Primer</b>	Forward: 5'-CGCAAATGGGCGGTAGGCGTG-3'
<b>Sequencing Primer</b>	Reverse: 5'-TAGTCAGCCATGGGGC GGAGA -3'
<b>Tags</b>	N/A
<b>Selection Marker</b>	Mammalian: Neomycin
<b>Format</b>	HSV Virus
<b>Appearance</b>	Liquid
<b>Storage Buffer</b>	Serum-free DMEM with 5% PEG6000
<b>Storage</b>	1 year when stored at -80°C.

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).