


## Lenti-UbC-GFP Control Lentivirus



**Cat. No.**  
LV027  
**Quantity**  
10ml

<b>Vector</b>	pLenti-III-UbC-GFP Control Vector
<b>Titer</b>	10 <sup>6</sup> IU/ml
<b>Selection Marker</b>	Mammalian: Puromycin
<b>Format</b>	Lentivirus
<b>Appearance</b>	Clear Liquid
<b>Storage</b>	1 year when stored at -80°C or lower in a non-frost free freezer.

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications.  
Please contact a technical service representative for more information (1-866-757-2414).